


# COVID-19 Ag Respi-Strip

## Detection test of SARS-CoV-2

1. In the test tube the sample is diluted 1:2 with the LY-S buffer
2. The total volume must be 200 $\mu$ l (buffer 100 $\mu$ l - sample 100 $\mu$ l)
3. Dip the strip (the liquid must not exceed the red line under the red arrows)
4. Leave to react for 15 minutes
5. Read the final result



### COVID-19 Respi-Strip Results

Any signal visible at 15 min on the test line, even a weak one, must be interpreted as positive

# COVID-19 Ag Respi-Strip

## Performances:

The kit was validated in comparison with RT-PCR on a total of 231 nasopharyngeal swab specimens in two Reference Hospitals: University Hospital Laboratory of Brussels (Evaluation 1) and the University Hospital Laboratory of Liège (Evaluation 2).

The following results were obtained:

<b>Evaluation 1 (n=99)</b>	(95% confidence interval)	
<b>Sensitivity</b>	60%	(40.7 to 76,8%)
<b>Sensitivity threshold</b>	85.7% on sampling with CT under 25 (62.6 to 96.2%)	
<b>Specificity</b>	100%	(93.4 to 100%)
<b>Positive predictive value</b>	100%	(78.1 to 100%)
<b>Negative predictive value</b>	85.2%	(75.2 to 91.8%)
<b>Agreement</b>	87.9%	(87/99)

<b>Evaluation 2 (n=132)</b>	(95% confidence interval)	
<b>Sensitivity</b>	60.3%	(48.1 to 71.3%)
<b>Sensitivity threshold</b>	76.7% on sampling with CT under 25 (61 to 87.7%)	
<b>Specificity</b>	98.3%	(89.7 to 99.9%)
<b>Positive predictive value</b>	97.8%	(86.8 to 99.9%)
<b>Negative predictive value</b>	66.7%	(55.7 to 76.2%)
<b>Agreement</b>	77.3%	(102/132)

## Detectability:

- Viral detectability :  $5 \times 10^3$  pfu/mL
- Recombinant protein detectability : 0,25 ng/mL

## Reactivity and cross reactivity:

<b>Viruses</b>	<b>Result</b>
Influenza A	neg
Influenza B	neg
Respiratory Syncytial Virus (RSV)	neg
Respiratory Adenovirus	neg
Parainfluenza	neg
Rhinovirus	neg
Metapneumovirus	neg
Enterovirus	neg
Coronavirus HKU1	neg
Coronavirus OC43	neg
Coronavirus 229E	neg
Coronavirus NL63	neg
Coronavirus SARS-CoV	pos
Coronavirus SARS-CoV-2	pos

<b>Bacteria</b>	<b>Result</b>
<i>Staphylococcus aureus</i>	neg*
<i>Legionella pneumophila</i>	neg
<i>Nocardia asteroides</i>	neg
<i>Streptococcus pneumoniae</i>	neg
<i>Moraxella catarrhalis</i>	neg
<i>Streptococcus pyogenes</i>	neg
<i>Haemophilus influenzae</i>	neg
<i>Pseudomonas aeruginosa</i>	neg
<i>Acinetobacter baumannii</i>	neg
<i>Klebsiella pneumoniae</i>	neg
<i>Mycoplasma pneumoniae</i>	neg
<b>Fungus</b>	<b>Result</b>
<i>Aspergillus niger</i>	neg

\* Staphylococcus aureus was found positive at high bacteria concentrations (10e9 cfu/mL). No false-positive results were obtained when testing naso-pharyngeal sample from known S. aureus infected patients.

## Coris BioConcept

Science Park CREALYS  
Rue Jean Sonet 4A  
5032 Gembloux - Belgium  
Ph: +32 (0)81 719 917  
Fax: +32 (0)81 719 919  
info@corisbio.com - sales@corisbio.com  
www.corisbio.com



D-CUS-10-23-S/03  
April 2020